Prodigy[®] No Coding Test Strip


NO CODING REQUIRED PRODIGY NO CODING TEST STRIP For Reference Numbers 52800, 52810, 52835, 72500, 73200, and 51800

- For in vitro diagnostic use only.
- All Prodigy® meters, Prodigy Lancing Devices, and Prodigy No Coding Test Strips" are for single patient use only.
- Do not share them with other people. Do not share them with family members!
- · Do not use on more than one patient!
- All parts of the kit should be treated as biohazardous and can transmit infectious diseases, even after you clean and disinfect.

Please read the owner's manual that came with your Prodigy® meter for steps to clean and disinfect.

Patient Cautions:

IVD * For in-vitro diagnostic use (external use only).

This system (meter) should not be used on critically ill nations.

PATIENT CAUTIONS

single patient use.

- 1. Do not use for neonates (newborns or infants).
- Do not use to screen or diagnose diabetes mellitus.
 Prodigy® Blood Glucose Monitoring Systems are for
- 4. Do not share them with other people!
- 5. Do not share with family members!
- 6. Do not use on more than one person!
- 7. Do not use for critically ill persons!
- 8. Do not use Alternate Site Testing (AST) to calibrate Continuous Glucose Monitors (CGM)!
- 9. Do not use AST to measure insulin doses
- 10. For in-vitro (external) diagnostic use only.
- Treat all parts of the kit as biohazardous that can transmit infectious diseases, even after you clean and disinfect.

NTENDED USE

- Your Prodigy meter is intended to be used for patient selfmonitoring and should not be used to collect blood from more than one person as this poses a risk of transmitting blood-borne pathogens such as Hepatitis B or HIV.
- The Prodigy® meter, Prodigy No Coding Test Strips and Prodigy® Lancing Device are for single patient use. Do not share them with anyone including other family members! Do not use on multiple patients!
- Treat all parts of the kit as biohazardous that can transmit infectious diseases, even after you clean and disinfect.
- Review the owner's manual that you received with your Prodigy® meter for steps to clean and disinfect.

All Prodigy® Blood Glucose Monitoring Systems are intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the finger tip, forearm, palm, upper arm, calf or thigh. All Prodigy® Blood Glucose Monitoring Systems are intended to be used by a single person and should not be shared.

All Prodigy® Blood Glucose Monitoring Systems are only for in vitro (outside the body) diagnostic self testing at home by people with diabetes. Self testing should be used as an aid to measure and review the effects of diabetes control. The Prodigy® Blood Glucose Monitoring System should not be used to diagnose or screen for diabetes, or for neonates (newborns or infants). AST should be done only during stable times when glucose is not changing quickly.

done only during stable times when glucose is not changing quickly. The Prodigy[®] No Coding Test Strips are for use only with a Prodigy[®] meter to measure glucose (sugar) in fresh capillary whole blood samples drawn from the finger tip, forearm, palm, upper arm, calf or thigh.

STORAGE AND HANDLING

To ensure that your Prodigy® No Coding Test Strips work as designed, please adhere to the sections that follow:

Important: Do not use Prodigy® No Coding Test Strips after the date

they expire or results may be wrong.

Prodigy® No Coding Test Strips expire 90 days after opened.

When you first open the vial, write the date opened on the

- When you first open the vial, write the date opened on the Prodigy® No Coding Test Strips vial.

 Store the Prodigy® No Coding Test Strips in a cool, dry place
- Store the Prodigy® No Coding Test Strips in a cool, dry place between 39°-104°F (4°-40°C) and between 10% and 85% relative humidity. Keep the Prodigy® No Coding Test Strips away from direct sunlight.
- Prodigy® No Coding Test Strips must be stored ONLY in

eir vial.

Do not transfer them to a new vial or other holder. Do not touch the Prodigy® No Coding Test Strips when your hands are wet.

- Use each Prodigy® No Coding Test Strips promptly after taking it from the vial. Close the vial lid quickly after taking a new strip.
- Keep the vial lid closed at all times when not in use.
- Do not bend, cut or alter the Prodigy® No Coding Test Strips.

HOW IT WORKS

Prodigy® No Coding Test Strips provide a measure of glucose in whole blood. The blood is drawn into the channel of the Prodigy® No Coding Test Strips that absorbs it through capillary action. The blood glucose level is measured and displayed on your Prodigy® meter. Only 0.7 μ L (microliters) of blood is required and the time to get results is 7 seconds.

FOR ACCURATE RESULTS

Prodigy® No Coding Test Strips should be used with a fresh blood sample. For best results read the sections below, along with your Prodigy® meter Owner's Manual before you obtain a blood sample. To get a drop of blood, follow these steps:

- Wash your hands and the puncture site. Rinse and dry yourself before testing. Do not clean the testing site with alcohol. Instead, wash the test site with soap and water only.
- Lance the puncture site to obtain a drop of blood. You may massage or gently rub the site in order to obtain the required amount of blood. Do not squeeze the puncture site in excess.

NO CODING REQUIRED PRODIGY® No Coding Test Strips

No calibration needed. Simply insert the Prodigy® No Coding Test Strips and proceed to the Control Solution Testing or Blood Glucose Testing sections of this insert.

CONTROL TESTING

Always perform control tests:

- · When you begin using a new vial of Prodigy® No Coding Test Strips.
- At least once a week
- If you suspect that your Prodigy® meter or Prodigy® No Coding Test Strips are not giving proper results.
- When your blood glucose test results are not consistent with how you feel, or when you think your results are not correct.
- If you drop the meter.
- When you change the batteries.

Prodigy® Control Solution is used to check testing technique and system performance. Prodigy® Blood Glucose Monitoring Systems use both high and low control solutions. When the control solution test result falls within the control range printed on the label of the vial of Prodigy® No Coding Test Strips being used, the system is working as it should, and your testing technique is good. Read your Prodigy® meter Owner's Manual for complete testing details.

Important: The control range can change with each new vial of Prodigy® No Coding Test Strips. Always use the control range on the label of your current vial of Prodigy® No Coding Test Strips.

BLOOD GLUCOSE TESTING

- Fully insert the Prodigy® No Coding Test Strips into your Prodigy®
 meter. Ensure that you insert the end of the strip with the contact
 bar into the test strip slot. The meter will turn on after you fully
 insert the Prodigy® No Coding Test Strips until it won't go further.
 The Prodigy® meter will perform several self-checks.
- 2. Follow the steps in your Prodigy® meter to obtain a blood sample of at least 0.7 µL of blood. Use a Prodigy® Lancing Device and Prodigy® Lancet to get the blood sample. Hold the blood drop to the end of the Prodigy® No Coding Test Strips that is outside the meter until blood has filled the confirmation window. The meter will been and begin to count down.
- After 7 seconds, your Prodigy® meter will display the blood glucose test result. The Prodigy® meter will store the test result. The meter will shut off when you remove the Prodigy® No Coding Test Strips. Discard the used Prodigy® No Coding Test Strips.

Test Results:

Your blood glucose test results are displayed either in milligrams per deciliter (mg/dL) or in millimoles per liter (mmol/L). The correct unit of measure in the USA is mg/dL. When your first purchase your system, please check to make sure the unit of measure is set to mg/dL. If it is not, please contact Customer Care at 1.800.243.2636.

Your Prodigy® meter is able to display test results from 20 to 600 mg/dL (or from 1.1 to 33.3 mmol/L). If your test result is below 20 mg/dL (1.1mmol/L), a warning message will appear showing a low glucose level. This may be a sign of severe low blood glucose. Treat this condition right away, as directed by your doctor or nurse.

Although this message could be due to a test error, it is safer to treat first and then do one or more tests.

High Glucose Values

If your Prodigy® No Coding Test Stripsresult is above 600 mg/dL, repeat your test with a new strip to confirm this reading. If the results should appear again, call your doctor or nurse right away.

Reference Values

Please consult your physician to establish your target ranges. The American Diabetes Association (ADA) suggests the following targets for most nonpregnant adults without diabetes. More or less stringent goals may be right for you.

Glycemic control Time of Day	ADA Suggested Reading
Preprandial plasma glucose (between meals)	70-100 mg/dL
Postprandial plasma glucose (1-2 hours after starting a meal)	<140 mg/dL

American Diabetes Association. Standards of medical care in diabetes – 2012. Diabetes Care. 2012;35 (Supp 1): S11-S63.

IF YOU SHOULD GET UNEXPECTED RESULTS

If your test results do not match how you feel:

- Make sure that the drop of blood fills the confirmation window of your Prodigy® No Coding Test Strips.
- Make sure that the Prodigy® No Coding Test Strips have not expired.
 Do not use Prodigy® No Coding Test Strips after the date they expire,
 which is listed on the label of the vial. Do not use Prodigy® No
 Coding Test Strips more than 90 days after you opened the vial of
 Prodigy® No Coding Test Strips.
- Conduct a Prodigy® Control Solution test to check that your Prodigy® meter and Prodigy® No Coding Test Strips are giving correct results.

High or low blood glucose levels may be a sign of a grave medical problem. If you get results that are high or low extremes, call your doctor or nurse.

Chemical components in sensor:

- Glucose Oxidase (Asperilligus niger) 10%
- Electron Shuttle 50%
- 3. Enzyme Protector 8%
- Non-reactive ingredients 32%
 TEST PRINCIPLE

A small electrical current is produced when glucose on the blood sample mixes with special chemicals on the Prodigy® No Coding Test Strips. This current is measured by your Prodigy® meter. Your Prodigy® meter converts the measure of current into a measure of blood glucose. The Prodigy® meter displays the measure of blood glucose.

Limitations:

- Hematocrit: Hematocrit levels below 20% or above 60% can cause false results. Please consult your doctor or nurse if you do not know your hematocrit level.
- Neonatal Úse: Prodigy® No Coding Test Strips should not be used for the testing of neonates (newborns or infants).
- Metabolites: Reducing substances such as uric acid and ascorbic acid when occurring in expected blood concentration, or acetaminophen, dopa, methyldopa, L-dopa and tolbutamide occurring in therapeutic blood concentration will not affect results.
- Lipemic Effects: Elevated blood triglycerides up to 2000 mg/dL do not significantly affect the results, however the blood glucose test may be affected beyond these levels.

INTERFERING SUBSTANCES

The Prodigy Preferred® Blood Glucose Monitoring System was tested for interfering substances. All the bias of test results were <10% compared with the measurements of the controlled pool. No obvious interference was observed in the interfering substance at therapeutic or physiological levels at three blood glucose levels. Severe dehydration and excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.

Elevated blood triglyceride, Reducing substances such as uric acid and ascorbic acid, Acetaminophen, Dopa, Methyldopa, L-dopa, and Tolbutamide (when occurring in normal blood or normal therapeutic concentrations) do not significantly affect results. However, abnormally high concentrations in blood may cause inaccurately high results.

Alternate Site Testing (AST)

- Use alternate site testing (AST) only two hours or more after a meal, two hours or more after taking insulin, or two hours or more after exercise.
- Do not use AST if you are aware that your glucose level is not as stable as usual.
- Do not use AST if you think you have hypoglycemia (low blood sugar).
- Do not use AST if you think you have hyperglycemia (high blood sugar).
- Do not use AST if your AST results do not match the way you feel.
- **6.** Do not use AST results to measure doses for insulin.
- Do not use AST results to calibrate continuous glucose monitors (CGMs).


Important: Using AST may lead to different results from each site. You may obtain a blood sample from your finger tip, palm, forearm, upper arm, calf or thigh. To obtain best results, and always if you suffer from hypoglycemia, you should test from the finger tip. Please read your Prodigy® meter Owner's Manual for more information about AST. Consult your doctor or nurse before you begin testing with a blood sample obtained from any site other than your finger tip.

PERFORMANCE CHARACTERISTICS Introduction to Linearity Study

each lot of test strips for a total of 390 tests.

To assess the linearity of Prodigy® Blood Glucose Monitoring Systems, a clinical study was performed using Prodigy® meters and three different lots of Prodigy® No Coding Test Strips. Different glucose concentrations were included in the study, ranging from 15 mg/dL. to 620 mg/dL. Readings were taken with the Prodigy® meter and the YSI 2300 blood glucose analyzer. Each Prodigy®

meter performed one test at each glucose concentration with

Regression analysis of the measurements using Prodigy® Blood Glucose Monitoring System and Glucose Analyzer YSI 2300.

Data Summary

Lot	Slope	Intercept	R2	r
All	0.9822	2.8123	0.9980	0.999
Lot 1	0.9757	4.1230	0.9980	0.999
Lot 2	0.9702	4.9569	0.9979	0.999
Lot 3	0.9810	3.1466	0.9980	0.999

Introduction to Alternate Site Accuracy Study To assess the accuracy of Prodigy® Blood Glucose Monitoring

70/89 (79%)

Systems on alternate sites, a clinical study was performed using Prodigy® meters, 3 different lots of Prodigy® No Coding Test Strips, and 100 participants. Different glucose concentrations were included in the study, ranging from 44 mg/dL to 449 mg/dL. Readings were taken with Prodigy® meters and the YSI 2300 blood glucose analyzer. Each Prodigy® meter performed one test from each site for each participant for a total of 600 tests The results for each site are presented below.


Fingertip test

38/89 (43%)

For glucos	e concentration	<75 mg/dL		
Within ±5 mg/dl (Within ±0.28 m			Vithin ±15 (Within ±0.8	
10/11 (91%)	11/11 (100%) 1	1/11 (100	%)
For glucose concentration ≧75 mg/dL				
Within +5%	Within +10%	Within +15%	. Withi	n +20%


89/89 (100%)


89/89 (100%)

Regression plot of fingertip capillary blood test 1 versus YSI 2300.

Palm test

Regression plot of palm capillary blood test 1 versus YSI 2300.


61/89 (69%)

Forearm test

33/89 (37%)

For gluco	se concentration	<75 mg/dL		
Within ±5 mg/c (Within ±0.28 n		±10 mg/dL ±0.56 mmol/L)		±15 mg/dL ±0.83 mmol/L)
5/11 (45%)	9/11 (8	2%)	11/11	(100%)
For glucose concentration ≧75 mg/dL				
Within ±5%	Within ±10%	Within ±15	5%	Within ±20%

82/89 (92%)

Regression plot of forearm capillary blood test 1 versus YSI 2300.

For alucose concentration <75 ma/dL

Upper arm test

 Within ±5 mg/dL (Within ±0.28 mmol/L)
 Within ±10 mg/dL (Within ±0.56 mmol/L)
 Within ±15 mg/dL (Within ±0.83 mmol/L)

 5/11 (45%)
 8/22 (73%)
 11/11 (100%)


For glucose concentration ≧75 mg/dL			
Within±5%	Within±10%	Within±15%	Within±20%
30/89 (34%)	56/89 (0%)	74/89 (83%)	86/89 (97%)

 $\label{eq:Regression} \textit{Regression plot of upper arm capillary blood test 1 versus YSI 2300.}$

For glucose concentration <75 mg/dL

Calf test

Within ±5 mg/dL (Within ±0.28 mr		Within ±1 (Within ±	0 mg/dL 0.56 mmol/L)		n ±15 mg/dL n ±0.83 mmol/L)
4 /11 (36%)		11/11 (1	00%)	11/11	(100%)
For glucos	e conc	entration	≧75 mg/dL		
Within ±5%	With	n ±10%	Within ±15	%	Within ±20%
28/89 (31%)	62/8	9 (70%)	79/89 (899	6)	89/89 (100%)


Regression plot of calf capillary blood test 1 versus YSI 2300.

Thiah test

87/89 (98%)

For glucose concentration <75 mg/dL Within ±5 mg/dL Within ±10 mg/dL Within ±15 mg/dL (Within ±0.28 mmol/L) (Within ±0.56 mmol/L) (Within ±0.83 mmol/L) 6/11 (55%) 10/11 (91%) 11/11 (100%) For glucose concentration ≧75 mg/dL Within ±20% Within +5% Within ±10% Within ±15% 33/89 (37%) 60/89 (67%) 74/89 (83%) 86/89 (97%)

Regression plot of thigh capillary blood test 1 versus YSI 2300.

Prodigy® Blood Glucose Monitoring System has been certified to meet the following standards: 98/79/EC,IEC 60601-1, IEC 61010-1, IEC 610

Please refer to your Prodigy® meter Owner's Manual for more detailed directions and illustrations.

PRODIGY*

Manufactured for: Prodigy Diabetes Care, LLC PO Box 481928

For Customer Care 24/7 call: 1.800.243.2636 www.prodigymeter.com

© 2014 Prodigy® is a registered trademark of Prodigy Diabetes Care.

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All rights reserved. TSMI Rev 8 02/14
Prodiay Manufacturado para Prodigy Diabetes Care LLC. 2701 Hutchison McDonald Road, Suite A Charlotte, N.C. 28269 .USA